# Tapered SwissPlus® Implant System

# Versatility Above And Below The Bone





# Overview

# The Tapered SwissPlus Implant System Simplifies Your Implant Cases By Providing A Comprehensive Solution For Extensive Restorative Applications.

The Tapered SwissPlus Implant is designed to be placed at tissue level in a one-stage surgical procedure. The fluted machined neck functions as the transmucosal extension of the implant receiving the prosthetic component of the restoration. The MTX® micro-textured titanium portion of the implant which includes the threaded area is placed subcrestally.

### The Tapered SwissPlus Implant System Has The Following Advantages:

- Only one surgical procedure is required, saving time for the clinician.
- The elimination of a second-stage surgery means less trauma for the patient.



# Prosthetic Components

## Impression Transfer, Cement-Retained Restorations\*

#### Impression-Taking

#### Indirect Transfer (w/ screw for Closed Tray Impressions)







Direct Transfer (optional procedure for Open Tray Impressions, Transfer supplied with Implant, screw sold separately)







**Implant Replica** 





## Cement-Retained Restorations

#### Straight Abutment/ Transfer (w/screw)

Narrow







Straight
Abutment
(w/transfer/waxing
coping) and Replica
(sold separately)







20° Angled Abutment for 6 positions (w/screw)







## Custom Restorations, Screw-Retained Restorations\*

Custom Restorations: "Cast-To" and **Plastic Castable** Abutments w/Screw

Single-Unit Restorations (Gold, Engaging)

Multi-Unit Restorations (Gold, Non-Engaging)

Multi-Unit Restorations (Plastic, Non-Engaging)

Ø 3.8 mm Platform SPGA w/MHLAS and OPS optimal, long waxing screw: MTWSD

N/A

N/A



OPGC w/GPCAS and OPS

OPCC w/GPCAS

Screw-Retained Restorations (for Ø 4.8 mm platform only)

**Tapered Abutment** Multi-Unit only (w/ healing cap)

**Tapered Abutment** Components Ø 4.8 mm Platform



OPACT w/TATHC





Indirect Transfer



ACTDT Tapered Abutment Direct Transfer



ACTR Tapered Abutment Replica Represents Tapered Abutment Attached to Implant



ACTGC Tapered Abutment Gold Coping w/Screw (SCTS) and Sheath (OPS)



ACTP 5.0 mm Plastic Castable Coping Includes Screw (SCTS)



TGC3 3.0 mm Tapered Abutment Bar Gold Copings w/

Screw (SCTS)



TGC5 5.0 mm Tapered Abutment Bar Gold Copings w/ Screw (SCTS)



TTC5 5.0 mm Tapered Abutment Titanium Bar Coping w/Screw (SCTS)



ACTT Titanium Temporary Includes Screw (SCTS)



## **Overdenture Restorative Components\***

#### Overdenture Restorations

Locator® **Abutment** 

## Ø 3.8 mm Implant Platform



Height



Height







4.0 mm Cuff Height Height

## Ø 4.8 mm Implant Platform



OPLOC4/1 OPLOC4/2 2.0 mm Cuff 1.29 mm Cuff Height



OPLOC4/3 3.0 mm Cuff Height



OPLOC4/4 4.0 mm Cuff Height

#### Locator Components

Note: Replacement Males are available. See product catalog.



LOCMP2 Male Processing Package



LOCIMP Impression Coping



LOCAN/4 Abutment Analog (for 3.8 mm Platform)



LOCAN/5 Abutment Analog (for 4.8 mm Platform)



LOCPP Parallel Post



LOCBF2 Bar Female



Cast-To Bar Female



LOCLBF Laser Bar Female

### Overdenture Restorations

Ball Abutments w/ **Cap Attachment** 



CAH & CAN included



OPBA3 CAH & CAN included



CAH Cap Attachment Housing




CAN Nylon Liner

### **Ball Transfer** Components



## SwissPlus Surgical Tools\*



## Prosthetic Tools and Torque Wrenches\*



 $<sup>^{*}\ \ \</sup>mathsf{Parts}\ \mathsf{shown}\ \mathsf{are}\ \mathsf{representational}\ \mathsf{in}\ \mathsf{nature}\ \mathsf{and}\ \mathsf{do}\ \mathsf{not}\ \mathsf{display}\ \mathsf{the}\ \mathsf{exact}\ \mathsf{geometry}\ \mathsf{of}\ \mathsf{the}\ \mathsf{abutment}.$ 



## Tapered SwissPlus Implant System\*

### **Tapered Implants**

(available lengths: 8.0 mm, 10.0 mm, 12.0 mm, 14.0 mm); includes 2mmL cover







Platforms with Fixture Mount







Surgical Cover Screws







<sup>\*</sup> Parts shown are representational in nature and do not display the exact geometry of the abutment.



Contact us at 1-800-342-5454 or visit

zimmerbiometdental.com

Zimmer Biomet Dental Global Headquarters 4555 Riverside Drive Palm Beach Gardens, FL 33410 Phone: +1-561-776-6700

of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authori Locator is a trademark of Zest IP Holdings, LLC. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0186 REV B 05/18 ©2018 Zimmer Biomet. All rights reserved.

